



SUGGESTED EXPOSURE FACTORS/TECHNIQUES

40 INCH SID	Small					AVERAGE			Large		
HAND	<i>(including children 5 yrs & up)</i>										
VIEWS	mA	Pulse Range			kVp*	mA	Pulses	kVp	mA	Pulses	kVp
PA/AP	15	8	-	12	54	15	16	56	15	16	62
OBLIQUE	15	8	-	12	56	15	16	58	15	16	64
LATERAL	15	8	-	12	56	15	16	60	15	16	66
WRIST											
VIEWS	mA	Pulse Range			kVp*	mA	Pulses	kVp	mA	Pulses	kVp
PA	15	8	-	12	54	15	16	58	15	16	64
OBLIQUE	15	8	-	12	56	15	16	60	15	16	66
LATERAL	15	8	-	12	56	15	16	62	15	16	68
ELBOW FOREARM	Small					AVERAGE			Large		
VIEWS	<i>(including children 5 yrs & up)</i>										
VIEWS	mA	Pulse Range			kVp*	mA	Pulses	kVp	mA	Pulses	kVp
AP	15	8	-	16	60	15	22	60	15	22	66
OBLIQUE	15	8	-	16	60	15	22	62	15	22	68
LATERAL	15	8	-	16	62	15	22	64	15	22	70

* May need to reduce kVp by 3 -5 on infants/children under 5 years of age