



The purpose of this document is to provide a step-by-step guide on how to create a new Send Destination or edit an existing Send Destination in the Opal PACS.

A Send Destination is established for sending study information to another Opal PACS/Non-Opal PACS or for off-site storage.


PLEASE NOTE: This functionality is a permission for those users in the ADMIN group only.

To begin:

1. Log into the Opal-Web PACS using your assigned user name/password combination.
2. Click on the  tab.



3. From the left menu, select **Send Destinations** from the DICOM category and the following screen will display:



You can perform the following functions:

1. **Edit** an existing destination
2. **Add** a new destination
3. **Delete** an existing destination

4. How to complete the fields:

TYPE

Select **OPAL** if sending to another Opal PACS

Select **CSTORE** if sending to any other non-Opal PACS

HOST NAME

Enter **IP Address** of Server you are sending to.

AE TITLE

Enter **Name of Box you are sending FROM.**

REMOTE AE TITLE

Enter **Name of Box you are sending TO**

PORT

*This is how the boxes talk and allows a communication portal

33336 if *OPAL to OPAL*

104 if *CSTORE*

DESCRIPTION

Enter send destination **name as you want it to appear** in drop down box when sending.

TRANSFER SYNTAX

Select resolution you want study to be sent in.

Little Endian

JPEG 2000

JPEG Lossless

START TIME


Indicate the time the particular study should be sent for this specific send destination queue

The screenshot shows a web interface for adding a new send destination. At the top right, there are two buttons: 'Add New' and 'Delete'. Below these is a form with the following fields and values:


- Type: OPAL (dropdown menu)
- Host Name: (empty text box)
- AE Title: (empty text box)
- Remote AE Title: (empty text box)
- Port: (empty text box)
- Description: (empty text box)
- Transfer Syntax: Little Endian (dropdown menu)
- Start Time: now (dropdown menu)

At the bottom of the form is an 'Add' button. Two blue arrows point from the text instructions on the left to the 'Type' and 'Description' fields in the form.

8. The following table will display with your list of created SEND DESTINATIONS.



9. Click **Send**



10. Click OK

11. The study will now be sent based on the criteria you indicated in that chosen Send Destination and will send at the specific time.