

DIGITAL X-RAY

MOBILE

SOLUTIONS

Revolutionary Imaging Systems



[2020imaging.net/mobile](http://2020imaging.net/mobile)

# Glassless DR Panel Choices

## MOMENTUM-DR GL-P FLATPANEL DR


2020imaging.net/momentum



- 14x17"
- Cesium Scintillator
- Lithium Ion Capacitor (7-yr warranty)
- Rapid Charge
- Wireless Interface
- Featherlight 3.9 lbs
- IP56 Dust & Liquid Resistance
- Easy Grip Design

## 14X17 UNIVERSAL DR FLATPANEL DR

Pending 510(k) Clearance



- 14x17"
- 16-bit A/D Conversion
- Cesium Scintillator
- Lithium Ion Battery (2 included)
- Wireless Interface
- Lightweight 6.6 lbs
- IP44 Dust & Liquid Resistance
- Carbon Fiber and Mg Alloy Construction

# Hospital-Grade Software



## ACQUISITION INTERFACE

- Under 6 second Processing
- Background Submission
- Preset X-ray Exams
- Adult & Pediatric Settings
- Unlimited Image Container
- Fully Customizable
- Advanced Image Processing
- Image Enhancement Tools
- Auto-Shutter & Image Enhance
- User Friendly Interface

### Multi-funton Studylist

- Auto DICOM Send Included
- Custom Search/Query
- Burn/Import CDs
- Customized User Settings



### Advanced Viewer

- Full Image Viewer on the go
- Fast, User Friendly Interface
- Image Fine Tuning
- Customized Layout



# Mobile SR-130 X-ray



## SR-130 FEATURES

- High Frequency Inverter
- Lower Radiation Dose\*
- Agile Lightweight Portability
- Hospital Quality Imaging
- Intuitive Interface
- Superb Image Quality
- Extended Tube Life
- UL & CUL Listed
- Built-in Cassette Storage

[2020imaging.net/sr130](http://2020imaging.net/sr130)

## SR-130 Specifications

- 3.0kw Output - Regulated mA/kVp
- KVP Range: 40 - 100
- MA Range: 15 & 30
- MAS Range: 0.15 - 120



\*Compared to standard frequency (single phase) x-ray.



Durable  
Superior  
Efficient  
Reliable  
Affordable

[2020imaging.net/mobile](https://2020imaging.net/mobile)

