



DIGITAL X-RAY

VETERINARY SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/vet



Wireless, Fixed, & Hybrid DR



MOMENTUM 14X17 WIRELESS FLATPANEL



- 3 Year Warranty
- 14"x17" Imaging Area
- Cesium Scintillator
- 100 µm | 5 lp/mm
- Wireless connectivity
- Cassette Sized
- 7-year lithium-ion capacitor



HYBRID 17X17 FLATPANEL



- 5 Year Warranty
- 17"x17" Imaging Area
- Cesium Scintillator
- 139 µm | 3.6 lp/mm
- Wireless and tethered connectivity
- Cassette Sized



WIRELESS 14X17 FLATPANEL

- 5 Year Warranty
- 14"x17" Imaging Area
- Cesium Scintillator
- 150 µm | 3.3 lp/mm
- Wireless connectivity
- Cassette Sized



TETHERED 17X17 FLATPANEL

- 5 Year Warranty
- 17"x17" Imaging Area
- Cesium Scintillator
- 139 µm | 3.6 lp/mm
- Fixed Tethered Interface
- Cassette Sized



Retrofit and New X-Ray



* RETROFIT ANY X-RAY

Wireless & Fixed DR Compatible



↳ VETSMART™ SYSTEM

- Fixed 40" SID
- Toshiba x-ray tube 140 HU 1-2mm focal spot
- kVP range: 40-150kV
- Single Phase 208-230VAC
- HF 32kW capacitor assisted 400mA generator
- 15" Touch Screen Operator Console
- Six exclusive table base/table top designs
- Push and open tray to access flatpanel
- Exclusive Push and Open Front panel latch mechanism for easy generator access

Veterinary Specific Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Submission
- Preset X-ray Exams
- Animal Specific
- Anatomical Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Adjustment Tools
- Touchscreen Enabled
- User Friendly Interface

Multi-funzion Studylist

- PM Software Integration*
- Burn/Import CDs
- Custom Search/Query
- Customized User Settings


Advanced Viewer

- Vertebral Heart Score Tool incl.
- Fast, User Friendly Interface
- Image Fine Tuning
- Exam Room Viewing



MOBILE VIEWING WITH M-OPAL¹

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search



CLOUD-BACKUP OPTIONAL

- Live Synchronization
- 448-bit Encryption (traffic)
- Enterprise Grade
- Fully Monitored
- Fully Customizable
- Recover Deleted Files



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/vet

