



DIGITAL X-RAY

CHIROPRACTIC

SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/chiro

Retrofit and New X-Ray

17X48 DR PANEL



- Wireless or Tethered Interface
- Cesium Scintillator
- Lithium Ion Battery
- 139 micron | 3.6 lp/mm
- Active Area: 17" x 48"

Our Universal-DR 17" x 48" single exposure, dose reducing, time saving full spine panel. Improving workflow efficiency while improving the patient experience.



DDR COMPATIBLE X-RAY THAT MOVES

SYSTEM SPECS DIFFER FROM STRAIGHT-ARM ABOVE

VERSATILE STRAIGHT-ARM

- Superior Image Quality
- 17x17 Fixed DR compatible
- Space-Saving Design
- Full Range of Motion
- Motorized Vertical Movement
- Detents on 30°, 60°, 90°, & 120°
- H.T. Cables
- 125Kvp, 300k HU, 0.6-1.2mm FS

2020imaging.net/straight-arm



Wireless, Fixed, & Hybrid DR

MOMENTUM 14X17 WIRELESS FLATPANEL



- Wireless Interface
- Rapid Charge
- Cesium Scintillator
- Lithium Ion Capacitor
- Featherlight 3.9 lbs
- Cassette Sized w/ Easy Grip
- Glassless DR Panel


HYBRID 17X17 FLATPANEL



- Wireless or Tethered Interface
- Cesium Scintillator
- Lithium Ion Battery
- Lightweight 9.5 lbs
- Cassette Sized
- 16-bit Images
- 139 micron | 3.6 lp/mm

WIRELESS 14X17 FLATPANEL

- Wireless Interface
- Cesium Scintillator
- Lithium Ion Battery
- Lightweight 8 lbs
- Cassette Sized
- 16-bit Images
- 150 micron | 3.3 lp/mm



TETHERED 17X17 DETECTOR

- Tethered Interface
- Cesium Scintillator
- 16-bit Images
- Cassette Sized
- 3.6 lp/mm




Symmetry-Sync Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Submission
- Preset X-ray Exams
- Adult & Pediatrics Available
- Anatomical Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Adjustment Tools
- Touchscreen Enabled
- User Friendly Interface

MULTI-FUNCTION STUDYLIST



- EHR Integration Ready*
- Burn/Import CDs
- Custom Search/Query
- Customized User Settings

ADVANCED CHIROPRACTIC TOOLSET

- Advanced Orthogonal Anal.
- Angles (Regular, Cobb)
- Gonstead Measurements
- Cervical Curve
- Blair
- Lumbar Curve
- George's Lines
- Orthospinology
- Ilium Analysis
- Full Spine Stitching
- Line Drawing
- Spine Labeling
- Measuring
- Penning Analysis


Examples of Full Spine Stitching & Gonstead Measurement ↗



ZERO FOOTPRINT WEB VIEWER

iPad®, iPhone®, Android™, tvOS, any mobile or smart device!

- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search
- Diagnostic DICOM Viewing



*Additional implementation fee may apply.

Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/chiro

