

DIGITAL X-RAY

MOBILE

SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/mobile

Wireless DR Panel Choices

Cesium Scintillator

Wireless Interface

Cassette Sized

14x17"

MOMENTUM WIRELESS FLATPANEL DR

2020imaging.net/momentum



16-bit A/D Conversion

Lithium Ion Capacitor
(7-yr warranty)

Featherlight 5.5 lbs

Rapid Charge

14X17 UNIVERSAL DR WIRELESS FLATPANEL DR



16-bit A/D Conversion

Lithium Ion Battery (2 Included)

Pixel Matrix - 2288 x 2784

Carbon Fiber and Mg Alloy Construction

3.3 lp/mm (CSI)

7.9 lbs (CSI)

Hospital-Grade Software



ACQUISITION INTERFACE

- Under **6 second** Processing
- Background Submission
- Preset X-ray Exams
- Adult & Pediatric Settings
- Unlimited Image Container
- Fully Customizable
- Advanced Image Processing
- Image Enhancement Tools
- Auto-Shutter & Image Enhance
- User Friendly Interface

Multi-funton Studylist

- Auto DICOM Send Included
- Custom Search/Query
- Burn/Import CDs
- Customized User Settings



Advanced Viewer










- Full Image Viewer on the go
- Fast, User Friendly Interface
- Image Fine Tuning
- Customized Layout



Mobile SR-130 X-ray



SR-130 FEATURES

-  High Frequency Inverter
-  Lower Radiation Dose*
-  Agile Lightweight Portability
-  Hospital Quality Imaging
-  Intuitive Interface
-  Superb Image Quality
-  Extended Tube Life
-  UL & CUL Listed
-  Built-in Cassette Storage

2020imaging.net/sr130

SR-130 Specifications

- 3.0kw Output - Regulated mA/kVp
- KVP Range: 40 - 100
- MA Range: 15 & 30
- MAS Range: 0.15 - 120



*Compared to standard frequency (single phase) x-ray.



Durable
Superior
Efficient
Reliable
Affordable

2020imaging.net/mobile

