



DIGITAL X-RAY

PODIATRY

SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/podiatry

Wireless and Tethered DR

NEW WIRELESS MOMENTUM DR PANEL

2020imaging.net/momentum



MFG BY
KONICA MINOLTA

TETHERED 1012V DR PANEL

2020imaging.net/1012v



📶 Wireless Interface	🔋 Lithium Ion Capacitor (7-yr warranty)
🍃 Featherlight 3.30 lbs	⚡ Rapid Charge
🔗 Cesium Scintillator	📊 16-bit Images
📏 100 micron 5 lp/mm	📏 10x12" Cassette Sized

🔌 Tethered Interface	🔌 Detachable & Durable Cable
📏 Lightweight 5.5 lbs	📏 10x12" Cassette Sized
🔗 Cesium Scintillator	📊 16-bit Images
📏 125 micron 4 lp/mm	

New and Retrofit X-Ray

NEW PXS710D X-RAY

- 🔊 High Frequency Output
- 🐱 Dose Efficient/Ext Tube Life
- 📺 Superb Image Quality
- ↔ Bi-directional Tube Movement
- 🔌 Integrated Cable Management
- 🛠 Enhanced Patient Support
- ⚠ UL Certified 450lbs weight limit
- 📄 Qualifies for ADA Tax Credit
- 🇺🇸 UL & CUL Listed
- 🇨🇦 Health Canada Approved

New X-Cel X-Ray models also available

RETROFIT ANY X-RAY

Wireless & Tethered DR Compatible
All options qualify for ADA Tax Credit



Example of older X-cel MB-700 modified with new 3.5" h SD Base and 1012V panel.

Hospital-Grade Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Background Submission
- Preset X-ray Exams
- Adult & Pediatric Settings
- Anatomical Program (APR)
- Fully Customizable
- Advanced Image Processing
- Image Enhancement Tools
- Touchscreen Enabled
- User Friendly Interface

ADVANCED OPAL-RAD VIEWER



Multi-functon Studylist ⌵ Advanced Viewer ➔

- EHR Integration Ready**
- Podiatric Toolset Available†
- Burn/Import CDs
- Fast, User Friendly Interface
- Custom Search/Query
- Image Fine Tuning
- Customized User Settings
- Exam Room Viewing



MOBILE VIEWING WITH M-OPAL†

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search



CLOUD-BACKUP INCLUDED**

- Live Synchronization
- 448-bit Encryption (traffic)
- Enterprise Grade
- Fully Monitored
- Fully Customizable
- Recover Deleted Files



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/podiatry

