



MOMENTUM PORTABLE
EQUINE
SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/equine



Momentum EQ All In One Portable Solution



Lightweight Design Under 18 lbs



Counterbalance for stability




Ergonomic Panel and Tablet Holders



Allows for easy handling



Custom Tablet Enclosure



All in one charging



Requires one power cable



Rugged Lightweight Aluminum Design



Scratch Resistant Finish

10x12 Premium DR Panel



Lightest panel & holder at 4lbs	Integrated Cable Management
Scintillator: Csl (Cesium)	Battery Life: > 150 images/4hrs
100 micron 4 lp/mm	Rapid Charge enabled
16-bit AD Conversion	Warranty: 3 years
Cassette sized	7 year lithium-ion capacitor

Veterinary Specific Software



• All-in-One Software Solution	• User-friendly Interface	• Clinical Review	• Automated Image Processing
• Equine Specific Protocols	• DR Console	• Image Storage	• Produces an optimized image

Advanced · Superior · Efficient · Reliable · Affordable



2020imaging.net/equine