




829 Virginia Rd. Ste A – Crystal Lake, IL 60014

7/29/2015

Please be advised you may be offered a free Windows 10 upgrade on your x-ray acquisition and server.

It is **not recommended to perform this update** as it **WILL** cause problems with the x-ray software, rendering it unusable! 2020 Support is in the process of testing Windows 10 with our software(s), and will provide an update our website once the update is safe to run on your system. For details on how to roll back to Windows 7/8.1 if you have upgraded to Windows 10, see ' [Rollback Windows 10](#).



Sincerely,

20/20 Imaging, a division of Konica Minolta Medical Imaging
866-734-6234 ext 1